

## APPLICATION OF GENOMIC TOOLS

ONE TECH TAKES IT ALL

- EVOLUTION IN TECHNIQUES
- TOOL BOX OF CGE GENOMIC EPIDEMIOLOGY

## SUB-TYPING TECHNOLOGIES 1920

SEROTYPING, PHAGE & PLASMID TYPING

PULSED-FIELD GEL ELECTROPHORESIS PFGE

MULTI-LOCUS SEQUENCE TYPING MLSI

NEXT-GENERATION SEQUENCING NGS

2 WEEKS -  
- 1 MONTH

≈ 6 DAYS

OLD	NEW
ISOLATION	ISOLATION DNA & LIB PREP NGS
CHARACTERIZATION & SUB-TYPING	BIOINFORMATICS PLUG & PLAY TOOLS

INTERPRETATION

ONE TECH THAT TAKES IT ALL SINCE  
IT WORKS FOR ALL KINDS OF BACTERIAS

SPECIES  
VIRULENCE  
PHYLOGENY  
RESISTANCE

} WHAT TO STUDY

## INTRODUCTION OF NGS


## NEXT GENERATION SEQUENCING

### LIBRARY PREPARATION

- FRAGMENT DNA
- LIGATE ADAPTORS ( CONTAIN AMPLIFICATION  
& SEQUENCING PRIMERS )

## • AMPLIFICATION

## • SEQUENCING



IF IN ONE RUN  
YOU HAVE MORE  
THAN 1 SAMPLE  
THEN CREATE A  
SEQUENCE  
SPECIFIC FOR EACH  
SAMPLE

## SEQUENCE DATA

STD DATA STORES → .FASTA FILES

- 1) HEADER + ID
- 2) SEQUENCE DATA (YOU CAN PUT MORE THAN ONE HERE)

>G1 | 297332897 | NEF ESCHERICHIA COLI  
1247 CHROMOSOME, COMPLETE

## RAW-READS

SEQUENCER OUTPUT IS NOT FASTA BUT FASTQ

FAST-Q = FASTA + QUALITY SCORES

- HEADER
  - DNA SEQ
  - NAME FIELD
  - QUALITY SCORES
- } 4 LINES PER READ

DE-MULTIPLEXING DONE BY SEQUENCERS

SINGLE-END READS

PAIR-END READ → FRAGMENT GAP

INSERT SIZE — ALL FRAGMENTS WITH NO ADAPTERS

## DATA QUALITY

- COVERAGE
- DEPTH
- BREADTH OF COVERAGE

TRIMMING DATA "N" NOT APPLICABLE  
REMOVING LOW-QUALITY DATA

DATA SPLITTING  
DATA CHIPPING  
DATA TRIMMING

## COVERAGE

$$C = N \frac{L}{G}$$

— READ LENGTH (points to L)  
— GENOME LENGTH (points to G)  
# OF READS (points to N)

## BREADTH OF COVERAGE

$$C_b = \frac{\text{ASSEMBLY SIZE}}{\text{TARGET SIZE}}$$

$$C_b \in [0, 1]$$

$P(\text{TO RECONSTRUCT WHOLE GENOME})$

### EXAMPLE

$$\begin{array}{l} N = 5 \cdot 10^6 \text{ READS} \\ L = 100 \text{ bp} \\ G = 5 \cdot 10^6 \text{ bp} \end{array} \left\{ \begin{array}{l} \text{COVERAGE} \\ C = 100 \times \end{array} \right. \begin{array}{l} 100 \text{ READS COVER} \\ \text{EACH POSITION} \\ \text{IN THE GENOME} \end{array}$$

## ASSEMBLY

RAW READS ARE VERY SHORT, SO WE MERGE THEM IN ORDER TO MAKE LONGER SEQUENCES

- **MAPPING TO REF.** (W/ REF. GENOME)  
STUDYING & ANALYZING GENE, SNP?

- 
- **DE-NOVO ASSEMBLY** (W/OUT REF. GENOME)  
DISCOVERING & IDENTIFYING NEW GENES
- 

ANALYSIS - MOST PROGRAMS IN UNIX  
IMPORTANT TO KNOW COMMAND-LINE TOOLS

NEW WEB BASE TOOLS  
EASIER TO USE

DIFFERENT SEQUENCES  
NEED DIFFERENT ASSEMBLERS

DE-BRUIJN GRAPHS

FROM READS TO CONTIGS  
(FASTQ) (PASTA)

### GOAL

IDENTIFY GENE OR  
GENETIC MARKER THAT  
YOU ARE LOOKING FOR

### GOOD CONTIGS

IF CONTIGS ARE LARGE  
ENOUGH TO CONTAIN A  
GENE OR COMPLETE MARKER  
THAT YOU ARE STUDYING

- **# OF CONTIGS** LONGER ARE BETTER  
EASY TO ASSEMBLE
  - **N50**
- ↓ YOU SHOULD WANT AS FEW AS POSSIBLE

MEDIAN SIZE OF YOUR CONTIGS

↑ THE BIGGER THE BETTER

SPLITTING BARCODE

TIME IS LOOK FOR QUALITY

ASSEMBLY (WITH OR WITHOUT REFERENCE)

DE NOVO ASSEMBLY

FROM RAW READS TO CONTIGS

TOOL DESCRIPTION & APPLICATIONS

N50 - HOW GOOD IS YOUR ASSEMBLY

MOST CGE ARE ASSEMBLY BASED